Hoefer TE 70 and TE 77

semi-dry transfer units

 **Hoefer**

Page finder

Semi-dry Transfer Unit function and description	1
Unpacking	2
Specifications	2
Important information.....	3
Operating instructions	4
Care and maintenance.....	10
Troubleshooting.....	11
Electrotransfer notes.....	13
Bibliography	15
Ordering information	16

Safety warnings and precautions

Important user information

Please read this entire manual to fully understand the safe and effective use of this product.

The exclamation mark within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance instructions in the literature accompanying the instrument.

The lightning symbol within an equilateral triangle is intended to alert the user to the risk of exposure to high voltages.

Hoefer, Inc. reserves the right to make changes in the specifications without prior notice.

Warranty and liability

Hoefer, Inc. guarantees that the product delivered has been thoroughly tested to ensure that it meets its published specifications. The warranty included in the conditions of delivery is valid only if the product has been installed and used according to the instructions supplied by Hoefer, Inc.

Hoefer, Inc. shall in no event be liable for incidental or consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related exposures, however caused, arising from the faulty and incorrect use of the product.

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Semi-dry Transfer Unit function and description

The TE 70 and TE 77 semi-dry transfer units rapidly transfer proteins from polyacrylamide gels onto a membrane by means of a low current and low voltage electrotransfer with minimal Joule heating. Most transfers are complete in one hour or less.

The smaller TE 70 unit transfer surface is 14×16 cm, and the larger TE 77 transfer surface is 21×26 cm. Several small gels of the same thickness can be laid side-by-side or two larger gels can be stacked into a carefully constructed multi-layered transfer stack. A plastic mask may be used to prevent current from bypassing the gel stack.


The base and cover body are made of inert Noryl® plastic. The base acts as a reservoir for buffer that may leak during the transfer. The anode in the base is platinum coated and the cathode in the cover is stainless steel. A safety interlock switch breaks the circuit when the cover is removed.

Fig 1. Semi-dry transfer unit main components.

A power supply with the following rating is required but not included:

TE 70: up to 200 mA, 50 V

TE 77: up to 450 mA, 50 V

This declaration of conformity is only valid for the instrument when it is:

- used in laboratory locations.
- used as delivered from Hoefer, Inc. except for alterations described in the user manual.
- connected to other CE labeled instruments or products recommended or approved by Hoefer, Inc.

Unpacking

Unwrap all packages carefully and compare contents with the packing list, making sure all items arrived. If any part is missing, contact Hoefer, Inc. Inspect all components for damage that may have occurred while the unit was in transit. If any part appears damaged, contact the carrier immediately. Be sure to keep all packing material for damage claims or to use should it become necessary to return the unit.

Specifications

TE 70	
Max. wattage	10 W
Max. voltage	50 V~
Max. amperage	200 mA
Gel size	up to 14 × 16 cm
TE 77	
Max. wattage	22.5 W
Max. voltage	50 V~
Max. amperage	450 mA
Gel size	up to 21 × 26 cm
Environmental operating conditions	Indoor use: 4–40 °C Humidity up to 80% Altitude up to: 2000 m Installation category: II Pollution degree: 2
Dimensions (w × d × h)	32 × 25 × 9 cm (12.6 × 9.9 × 3.6 in.)
Product certifications	EN61010–1, UL3101–1, CSA C22.2 1010.1, CE

English

Important information

- The electric components in the safety interlock housing of the transfer unit base must not become wet. Do not immerse any part of the unit in water. Rinse only the electrodes with distilled water before use. (Refer to the Care and maintenance section for cleaning instructions.)
- Be sure to apply enough buffer-soaked sheets of blotting paper under the membrane and over the gel to prevent the stack from drying out. Irreparable damage to the unit will result if the stack dries out and the unit is allowed to overheat.
- If this equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.

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Note: See Electrotransfer notes on page 11 for transfer tips as well as two buffer recipes and references.

Note: Always wear gloves when handling membranes to avoid leaving fingerprints.

Operating instructions

To transfer proteins, prepare the unit and transfer stack components, assemble the stack, and connect the leads to a power supply. Then run the transfer for the required amount of time. Each step is described below.

Assemble the transfer stack

1

Prepare the unit by rinsing the anode and cathode with distilled water, taking great care to keep water away from the safety interlock housing and leads.

2

Prepare the gel

Cut away the wells and/or stacking gel section. Equilibrate the gel in transfer buffer if required. Measure the gel.

3

Cut a mask to the proper size

Cut an opening in the solid mask to a size slightly smaller than the gel.

The required opening is centered and 2 mm smaller than the gel on all sides. Because the surface is slick and hard, take care when cutting mask with a tool such as a mat knife.

Place the mask in the base of the unit, centering the opening.

4

Prepare the blotting paper

For each gel, cut at least 6 pieces of blotting paper the *same size as the gel or slightly smaller*.

Gauge the thickness or number of blotting paper layers according to the amount of buffer required; up to 250 ml of buffer may be required for larger gels or for transfers of 60 minutes in order to prevent the stack from drying out or the buffer from becoming depleted. Saturate at least 3 pieces of the blotting paper with transfer buffer. One by one, center each sheet over the cutout on the mask and remove all trapped air by rolling a clean pipet from the center toward the edges. Check that the blotting paper completely covers the mask cutout and slightly overlaps the mask edges on all sides.

5

Prepare the membrane

For each gel, cut 1 membrane the *same size* as the gel or slightly smaller. (A larger membrane may contact an electrode panel, creating a pathway by which current can bypass the gel.)

Pre-wet nitrocellulose or nylon membranes with distilled water. Pre-wet PVDF or other hydrophobic membranes with methanol. Then soak all membrane types in transfer buffer for 2–5 minutes.


Important! Stack each layer with care, with edges parallel. As each layer is added, remove all air pockets by rolling a clean pipet from the center to the edges. Add a few drops of buffer to trouble areas to help remove air pockets.

Complete the stack

- a) Place the pre-wet membrane onto the stack of blotting paper.
- b) Place the gel on the membrane. Note: Proteins bind to the membrane as soon as contact occurs, so it is important to place the gel correctly on the first try.
- c) Cover the gel with several layers of saturated blotting paper.

Multiple gels: Either lay gels of the same thickness side-by-side (the position of the gels must coincide with cutouts in the mask), or stack 2 sandwiches layered as shown on the next page.

Fig 2. Transfer stack for a single gel.


Note: Take care to place the gel correctly on the first try because proteins begin to transfer immediately; once transfer has begun, moving the gel will distort results or cause “shadow” bands on the blot.

Fig 3. Transfer stacks for gels placed side-by-side.

Check electrical continuity within each stack: several sheets of buffer-soaked blotting paper must cover each gel, and each stack must fully contact both electrodes.

Note: There must be no electrical contact between the two stacks.


If two gels are stacked, separate them with porous cellophane—**not plastic wrap!** (Cellophane permits electric current to pass but stops proteins.) Cut the cellophane slightly smaller than the gel and wet with transfer buffer. Several sheets of buffer-soaked blotting paper on *each* gel provides electrical continuity.

Fig 4. Transfer stack for stacked gels.

Separate the stack for each gel with a sheet of cellophane.

Note: When transferring multiple gels, transfer efficiency depends on such factors as gel thickness, gel position in the stack, transfer buffer, membrane type, and, most importantly, the characteristics of the protein. The gel closest to the anode generally transfers the most completely, so it is preferable to lay gels side-by-side rather than stacking them.

Note: Staining the gel(s) for residual protein gives an indication of the completeness of transfer.

Note: Rewet dried membranes in the appropriate wetting buffer prior to processing.

7

Slide the 3 guides on the cover into the 3 grooves in the transfer unit base. Do not remove the cover until after the transfer is complete in order to prevent stack components from sliding out of alignment.

8

If transferring more than one gel, place a 1 kg weight on the cover to improve the electrical continuity of the stack. Note: A weight of more than 1 kg may press too much buffer out of the blotting paper, which will adversely affect the transfer efficiency.

Transfer

1

Set the power supply high voltage limit to 50 V

If possible, set the voltage limit before connecting the transfer unit. Use only an approved power supply such as the Hoefer PS2A200.

2

The power supply should be switched off and both the current and voltage controls turned to zero

Plug the short safety-interlock lead from the cover into the jack on the base. Then plug the color-coded leads from the base of the transfer unit into the power supply jacks, matching red to red and black to black. Do not reverse polarity. (Red is positive [+], and black is negative [-].)

3

Set the power supply current

The maximum current setting must not exceed 0.8 mA/cm^2 of the gel surface. If transferring several layers of gels, a longer transfer time may be required.

Use the graph to quickly find the maximum current setting for your gel size, or calculate the gel area (cm^2) and multiply it by 0.8 mA/cm^2 .

If available, set the power supply timer


Most transfers are complete within one hour, but larger proteins, proteins from native gels, and thicker gels may require an additional 1 hour of transfer time. The optimum transfer time for each protein and gel system must be determined empirically.

When developing a new protocol, it is important to monitor the voltage closely. It should not exceed 50 V. High initial voltage (20 V or more) indicates inappropriate ionic strength; check the buffer concentration and quality. If the voltage increases more than 10 V during the run, the stack may be drying out and scorching may result, which can damage the sample or the transfer unit. Another indicator of buffer depletion is increasing temperature; feel the top of the cover for excessive heat periodically, especially after one hour. If drying is indicated by either high voltage or increased temperature, stop the transfer and add more layers of buffer-saturated blotting paper.

Fig 5. Recommended maximum current settings for different gel sizes.

Note: Generally, smaller fragments transfer more quickly than larger ones.

Note: After 1 hour of transfer time, switch off the power supply and re-wet the stack. Transfer can then proceed for an additional hour.

After the transfer is complete

1

Turn off the power supply and disconnect the leads from the power supply.

2

Unplug the safety-interlock lead between the cover and the base of the unit. Remove the cover slowly because the stack may adhere to it.

3

Remove and dispose of the blotter paper.

4

Remove the gel(s). *Optional:* Stain for residual protein.

5

If desired, label the gel contact side of the membrane with a soft pencil.

6

Remove the membrane(s) from the stack with blunt forceps. Process the membrane following instructions accompanying your protocol or allow the membrane to air dry prior to storage.

7

Remove the remaining blotter papers and dispose. Remove the mask. Rinse the unit according to the Care and maintenance instructions.

Care and maintenance

- Do not autoclave or wash the unit in a dishwasher.
- Do not wet the safety-interlock housing or leads.
- Do not immerse the unit in water!

Rinse the cathode in the cover and the anode in the base with distilled water. Let the unit air dry completely. If using radioactive reagents, decontaminate the unit with a cleaning agent such as Decon™ 90 or Contrad™ 70. Never use abrasive cleansers.

Troubleshooting

problem	solution
Incomplete transfer	
<i>Blank or faint areas on the membrane</i>	<p>Remove trapped air pockets between the gel and membrane during stack assembly.</p> <p>Reduce the current setting to minimize the amount of gas produced by the electrolysis of buffer.</p> <p>Use buffer with a lower ionic strength.</p>
<i>Molecules do not migrate out of gel</i>	<p>Check all electrical connections. Check that current is flowing by plugging a multi-meter into the circuit and comparing the meter reading with the power supply reading.</p> <p>Check that the buffer pH is close to the intended pH. Most buffers should not be titrated; make fresh buffer.</p> <p>Use 3.5 mM SDS (0.1%) in the transfer buffer.</p> <p>Add several more sheets of buffer-saturated blotting paper to each side of the gel sandwich so that more buffer is present during the transfer.</p> <p>Do not use staining or fixing agents on the gel before transfer.</p> <p>Use a thinner gel.</p> <p>Reduce the gel acrylamide concentration.</p> <p>Increase the transfer period. Large fragments may require an additional hour.</p> <p>If using a non-nitrocellulose membrane, avoid including methanol in the transfer buffer or reduce the amount to the minimum possible.</p> <p>Use reagent-grade chemicals.</p> <p>Increase the net charge on the protein by using a transfer buffer with a different pH. Lower pH (<6–7) increases the positive charge on proteins; higher pH (>6–7) increases the negative charge on proteins.</p>

Smeared or diffuse band patterns	<p>If equilibrating before the transfer, shorten or eliminate the equilibration time and/or equilibrate under coldroom conditions.</p> <p>If the transfer buffer contains $\geq 10\%$ methanol, equilibrate the gel in transfer buffer for 30 minutes to allow it to shrink before assembling the stack. Note: Large proteins may not migrate as readily once the pore size is slightly reduced.</p> <p>Take care that the gel does not shift once it contacts the membrane.</p> <p>If excess heating occurs during the transfer, note if the stack is drying out. On subsequent transfers, either shorten the transfer time, reduce the current, or increase the number of buffer-soaked blotting papers in the stack.</p> <p>Check that the preferred binding surface of the membrane (if any) contacts the gel.</p>
Uneven band transfer	<p>The blotting paper and membrane must be the same size as the gel or 1–2 mm smaller. Larger sizes will provide an electrical path for current to bypass the gel.</p>
Inefficient binding to membrane	
<i>Chemical parameters</i>	<p>Fix or crosslink the molecule onto the membrane according to the requirements of the nucleic acid, protein, or membrane type.</p> <p>Prepare transfer buffer without SDS. (SDS can improve transfer efficiency but can interfere with protein binding to a PVDF membrane.)</p> <p>Add 10–20% methanol to the transfer buffer to enhance binding to nitrocellulose.</p>
<i>Membrane parameters</i>	<p>Use a membrane with a smaller pore size (0.10–0.20 μm) if proteins pass through the membrane.</p> <p>Place a membrane both over and under the gel to capture any proteins migrating in the opposite direction.</p> <p>Check if too much sample is available for the binding surface area by placing two membranes instead of one. If “blow through” occurs, reduce the sample load.</p> <p>Wear gloves when handling membranes.</p> <p>Store membranes at ambient temperature and out of direct sunlight.</p>

Electrotransfer notes

- Run the transfer as soon as possible after electrophoresis to minimize protein diffusion within the gel.
- Stacked gels must all be the same size.
- Limit transfers to two hours or less, rewetting the stack after 1 hour.
- The recommended methanol concentration for different membrane types are:

membrane type	methanol %
Charged nylon	0
Nitrocellulose	20
PVDF	15

- Use a buffer with low ionic strength such as one of the two listed below to prevent overheating. Use the CAPS buffer when Tris cannot be used (e.g. peptide sequencing). CAPS can improve transfer because of its effect on the charge of the protein (see Matsudaira, 1987).

Towbin buffer

(25 mM Tris, 192 mM glycine, 20% v/v methanol, pH 8.3, 1 liter)

Tris (FW 121.1)	25 mM	3.0 g
Glycine (FW 75.07)	192 mM	14.4 g
SDS* (FW 288.4)	0.1% (3.5 mM)	1.0 g

Dissolve in 600 ml distilled water. Add methanol as required[†]. Bring to 1 liter with distilled water.

Do not adjust the pH, which should be between 8.2 and 8.4. *Optional*: Chill before use.

**Optional*: Adding SDS can improve transfer efficiency.

[†]Depending on the membrane type selected (see table above), adding methanol can improve transfer results.

Note: Buffers containing methanol may deteriorate if stored for long periods—add methanol just prior to transfer.

CAPS buffer, 1X

(10 mM CAPS, pH 11.0, 1 liter)

CAPS (FW 221.3) 10 mM 2.2 g

[3-(cyclohexylamino)-1-propanesulfonic acid]

Dissolve in 600 ml distilled water, adjust to pH 11.0 with conc. NaOH. Adjust volume to 1.0 liter.

- For a 3-buffer system, refer to Kyhse-Anderson, J. (1984).
- Transfer efficiency varies depending on the gel concentration, which can be optimized. For more information, refer to Smejkal and Gallagher (1994).

Bibliography

- Bjerrum, O.J., Larsen, K., and Heegaard, N., *CRC Handbook of Immunoblotting of Proteins 1*, Section 7. CRC Press (1988).
- Gallagher, S., Winston, S.E., Fuller, S.A. and Hurrell, J.G.R., Immunoblotting and Immuno-detection. In *Current Protocols in Molecular Biology*. 10.8.1–10.8.17. Greene Publishing and Wiley-Interscience, NY (1993).
- Hancock, K. and Tsang, V., India ink staining of proteins on nitrocellulose paper. *Anal. Biochem.* **133**, 157–162 (1983).
- Kyhse-Anderson, J., Electroblothing of multiple gels: A simple apparatus without buffer tank for rapid transfer of proteins from polyacrylamide to nitrocellulose. *J. Biochem. and Biophys. Meth.* **10**, 203–209 (1984).
- Matsudaira, P., Sequence from Picomole Quantities of Proteins Electroblooded onto Polyvinylidene Difluoride Membranes. *J. Biol. Chem.* **262**, 10035 (1987).
- Sasse, J. and Gallagher, S., Detection of Proteins on Blot Transfer Membranes. In *Current Protocols in Molecular Biology*. 10.7.1–10.7.3. Greene Publishing and Wiley-Interscience, NY (1991).
- Smejkal, G., and Gallagher, S., Determination of Semidry Protein Transfer Efficiency with Transverse Gradient Gel Electrophoresis. *Biotechniques*. **16**, 196–202 (1994).
- Tovey, E. and Baldo, B., Comparison of semi-dry and conventional tank-buffer electrotransfer of proteins from polyacrylamide gels to nitrocellulose membranes. *Electrophoresis* **8**, 384–387 (1987).

Ordering information

product	quantity	code number
TE 70 Semi-dry Transfer Unit, 14 × 16 cm Includes 25 sheets of blotter paper, 50 sheets of cellophane, 2 uncut masks, precut mylar mask for 7 × 8 cm gels and for 14 × 16 cm gels.	1	TE70
Mask and Cellophane for the TE 70		
Porous cellophane, 20 × 35.5 cm	50	TE73
Uncut mask, 16.5 × 18.5 cm	4	TE74
TE 77 Semi-dry Transfer Unit, 21 × 26 cm Includes 25 sheets of blotter paper, uncut mask and cellophane. (Gels 20 × 25 cm do not require a precut mask.)	1	TE77
Mask and Cellophane for the TE 77		
Porous cellophane, 35 × 44 cm	50	SE1142
Uncut mask, 23 × 28 cm	4	TE78
Blotter Paper		
Blotter paper, precut, 14 × 16 cm	25	TE76-1146
Blotter paper, precut, 20 × 25 cm, for TE 77	25	TE76
Power Supply		
EPS 2A200	1	PS2A200

Noryl is a registered trademark
of General Electric Company.

Conrad 70, Decon and Neutrad are
trademarks of Decon Laboratories, Inc.

Printed in the USA

Hoefer, Inc.
953 Indiana Street
San Francisco, CA
94107 USA

www.hoeferinc.com